



## PACKAGE

Liquid / 1L (21H701A)

Liquid / 60mL (21H701D)

Liquid / 60mL x 2 (21H701D2)

Liquid / 60mL x 4 (21H701D4)

Liquid / 60mL x 6 (21H701D6)

ⓘ 2 - 8°C

IVD *In vitro* Diagnostics

⊗ **DO NOT** freeze

⊗ 18 months/block from light

## 1. PURPOSE OF USE

HbA1c Diluent is a hemolysate solution use for dilution whole blood sample or red blood cell for quantitative determination of Measure HbA1c II by Immunoassay.

For *In Vitro* Diagnostics use only.

## 2. STRUCTURE AND COMPOSITION OF KIT

A mixed solution of surfactants and some other additives.

- Good buffer pH 7.6 < 0.0001M
- Sodium Azide < 0.1%
- Water fill up to 100%

Bottle 1L / 60mL / 60mL x 2 / 60mL x 4 / 60mL x 6.

## 3. DIRECTION FOR USE

### *Preparation method of HbA1c Diluent.*

HbA1c Diluent is ready for use. Use it with Measure HbA1c II Kit. Please read insert carefully.

### *Storage*

It can be used until it expire under storage in refrigerator (2 - 8°C).

## 4. STORAGE AND SHELF LIFE

Storage : Store in cool and dark place (2 - 8°C).

Shelf Life: 18 months from manufactured date.

Expiry date is printed on the label of bottle.

## 5. PRECAUTION FOR USE AND HANDLING

- This product and it's components have toxicity. Take same precautions against infection as patient samples such as to wear gloves. Do not pipette by mouth.

- All samples and diluent, as well as their containers must be treated according to the handling manual for discarding infectious waste materials, or treat them according to the following methods.

• Sterilize with an autoclave at 121°C for 20 minutes, however do not process wastes containing sodium hypochlorite with an autoclave.

• Immerse at least one hour in sodium hypochlorite solution ( active chloride concentration of over 1000 ppm).

## 6. REFERENCES

1. M.J. York, in A Comprehensive Guide to Toxicology in Nonclinical Drug Development (Second Edition), 2017
2. In house data, UMA Diagnostics

## 7. MANUFACTURER

UMA Co., Ltd.

2-19-6 Yokosuka, Matsudo City, Chiba

Prefecture 270-0031

TEL: 047-710-4871 (dial-in)

FAX: 047-710-4872