

UMA CO., LTD.

2-19-6 Yokosuka

Matsudo, Chiba, Japan



MEASURE HDL A

Reagent for determination of HDL Cholesterol

Inhibition/Direct Method

↓ 2 - 8°C

IVD *In vitro* Diagnostics

QUALITY MANAGEMENT SYSTEM (BY TUV)

⊛ **DO NOT** freeze

⌚ 18 months/block from light

ISO 13485:2016

1. PURPOSE OF USE

Providing a quantitative *in vitro* assay for the High density lipoprotein cholesterol (HDL-c) concentration in serum or plasma.

2. GENERAL INSTRUCTION

- For *in vitro* diagnostics use only.
- Diagnosis should be made in a comprehensive manner, in accordance with other related test results and clinical symptoms by the doctor in attendance.
- For guaranteed results, usage of this product must comply with the instruction in this manual.
- If you use automatic analyzers, follow their instructions carefully.

SUMMARY

The five different types of lipoproteins are chylomicrons, very-low-density lipoprotein (VLDL), intermediate-density lipoprotein (IDL), low-density lipoprotein (LDL), and HDL. Lipoproteins classify according to their density and composition. Lipoproteins are complex particles that transport lipids, such as phospholipids, triglycerides, and cholesterol, between cells. HDL, as per its denomination, is the highest density of the lipoproteins, with the highest proportion of proteins to lipids. HDL is of particular interest in medicine, as research has shown a strong inverse association between HDL cholesterol concentration and the risk of atherosclerosis.

HDL is composed of cholesterol, triglycerides, and various apolipoproteins. In particular, the composition of HDL is apolipoproteins Apo-AI, Apo-AII, Apo-AIV, Apo-AV, Apo-CI, Apo-CII, Apo-CIII, and Apo-E.

The primary function of HDL is the transport of cholesterol from the peripheral tissues to the liver, playing a role in the biodistribution of lipids. HDL is known for its anti-atherogenic and anti-inflammatory properties, thanks to its uptake and return of the cholesterol stored in the foam cells of atherosclerotic plaques to the liver. Thus, reducing the size of the plaque and its associated inflammation.

3. MATERIALS REQUIRED BUT NOT INCLUDED

- Saline 0.9 % and high grade purified water
- Micropipet and other basic laboratory equipment.
- Lipids Calibrator and Lipids Control / MEASURE Human Lyo L-1 and MEASURE Human Lyo L-2

4. REAGENT COMPOSITION & PREPARATION

- Reagent R-1: ADPS; ascorbate oxidase, H₃PO₄.

Reagent R-1 is ready for use

- Reagent R-2: CHOD, Cholesterol Esterase (CHER), Peroxidase (POD), 4-aminoantipyrine (4-A-A), sodium azide

Reagent R-2 is ready for use

- Once open, Reagent stored on board the instrument is stable for 30 days with Hitachi 7180 Analyzers.
- Applicable to various automated analyzers.
- Calibrator Lipids Calibrator (separately sold): Put 2 mL of purified water to the vials of calibrator (Lipids Calibrator); leave at room temperature for 120 minutes and sometimes gently invert the vial before use. After reconstituting, Calibrator can be used without dilution.
- Control Lipids Control (separately sold): Put 2 mL of purified water to the vials of controls (Lipids Control); leave at room temperature for 120 minutes and sometimes gently invert the vial before use. After reconstituting, controls can be used without dilution.
- Controls MEASURE Human Lyo L-1 and MEASURE Human Lyo L-2 (separately sold): Put 5 mL of purified water to the vials of controls (Lyo L-1 and Lyo L-2); leave at room temperature for 45 minutes and sometimes gently invert the vial before use. After reconstituting, controls can be used without dilution.

c. Performance

- Sensitivity: Change in absorbance when using purified water is 0.001 - 0.050, and change in absorbance using standard solution (HDL cholesterol 1.295 mmol/L) as specimen is 0.05 - 0.17.

- Accuracy: When measuring a control sample, the result is within $\pm 10\%$ of assigned value.

d. Precision (on Biolis 30i / SK300)

Representative performance data on the analyzers are given below.

Results obtained in individual laboratories may differ.

Precision was determined using controls followed the CLSI Approved Guideline EP5-A2 with repeatability, reproducibility and total precision (1 aliquot per run, 2 run per day, 20 days). The following results were obtained.

Criterion: CV of Repeatability (aka. Within-run precision) is less than 3% and Total Precision is less than 5%.

Repeatability	Mean mmol/L	SD mmol/L	CV %
Control Lyo L-1	1.44	0.01	1.01
Control Lyo L-2	1.17	0.02	1.17

Reproducibility	Mean mmol/L	SD mmol/L	CV %
Control Lyo L-1	1.44	0.03	2.21
Control Lyo L-2	1.17	0.03	2.34

Total precision	Mean mmol/L	SD mmol/L	CV %
Control Lyo L-1	1.44	0.03	2.33
Control Lyo L-2	1.17	0.03	2.65

e. Correlation Test

Same measuring principle - Serum:

Regression equation: $y = 1.002x - 0.9$ ($n = 90$)

Correlation coefficient: $r = 0.999$

(y: Value obtained from using this method)

Reference Materials for Calibration

ReCCS JCCRM 223

10. EXPECTED VALUES

Normal reference range

- Male 0.98 - 2.33 mmol/L
- Female 1.24 - 2.75 mmol/L

Reference range should be established at each facility and judgement should be based on measurement results in a comprehensive manner together with clinical symptoms and other measurement results.

11. INTERFERENCES

- Icterus: No significant interference of conjugated/free bilirubin concentration up to 20 mg/dL
- Hemolysis: No significant interference of hemoglobin concentration up to 500 mg/dL
- Lipemia (Intralipid): No significant interference triglycerides concentration up to 3000 FTU
- Ascorbic Acid: No significant interference of ascorbic acid concentration up to 50 mg/dL
- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings. Please use another methods if the result is affected by any factors

12. HANDLING, USAGE & DISPOSAL**Handling**

1. Specimen can be potentially positive for infectious agents including hepatitis B virus and HIV. Wear glove and goggle when needed.
2. In case reagents got into skin, eye or mouth by mistake, wash it immediately with plenty of water and consult the doctor if needed.
3. If reagents are spilled, dilute with water and wipe it out. If specimen is spilled, spray 80% of alcohol over the specimen and wipe it out.

Usage

1. Store reagents under specified condition. Do not use after expiration date.
2. Do not use the container and auxiliaries included in this kit for other purposes.
3. Do not mix reagents of different lot for use.
4. Do not add to the reagent being used even if it is the same lot number.

Disposal

1. All specimens, as well as all instruments (e.g. test tubes) that come in contact with the specimens, must be treated by the following methods, or they must be treated according to the manual for infectious medical waste provided in each facility.

- Sterilize with an autoclave, subjecting them to high pressure saturated steam at 121 °C for more than 20 minutes. Do not process waste containing sodium hypochlorite solution with an autoclave.

• Immerse at least one hour in sodium hypochlorite solution (active chloride concentration of over 1000 ppm).

2. This reagent contains sodium azide. Sodium azide can react with lead pipe and/or steel pipe and can generate explosive metal azide. Make sure to use plenty of water at disposal. Concentration of sodium azide in R-2 is 0.05%.

13. INFORMATION FOR AUTOANALYZERS

❖ For Hitachi Model

Calculation Method		Two point
Temperature		37°C
	Specimen	2.0
Volume (μL)	R1	150
	R2	50
Wavelength (nm)	Main	546
	Sub	800
Measurement (cycle)	Point 1	10
	Point 2	16
	Point 3	34
Calibration type		Linear
Unit		mg/dL

14. OTHER INSTRUCTIONS AND CAUTION

- Results may differ depending on the sample/reagent ratio. Adjust parameters for different analyzer.

- Perform the QC procedure on the day of determination.

15. PACKING AND KIT CONFIGURATION

Code	Package	Test/Kit*	Test/Kit**
11H022B	1x30mL; 1x10mL	160	250
11H022A	1x60mL; 1x20mL	310	540
11H022A2	2x60mL; 2x20mL	620	1080
11H022A3	3x60mL; 3x20mL	930	1620
11H022A4	4x60mL; 4x20mL	1240	2160
11H032A	5x60mL; 5x20mL	1550	2700
11H022A6	6x60mL; 6x20mL	1860	3240
11H022	1x90mL; 1x30mL	470	810
11H022-2	2x90mL; 2x30mL	940	1620
11H032	3x90mL; 3x30mL	1410	2430
11H022-4	4x90mL; 4x30mL	1880	3240
11H022-5	5x90mL; 5x30mL	2350	4050

* For middle-scale automatic analyzers such as: SK300; BS series; BA200; BA400, Chemwell Series; Dirui Series; Biolzyer series, HumanStar 300, Erba Series; Bioelab Series, BX 3010; Pictus P500;...

** For large-scale automatic analyzers such as: CA800; CA400; Randox Imola; Randox Modena+; BM 6010; Biolis50j; SK500; AU Series; Pictus P700; C series; C1 series; HumanStar 600; Kenolab series ...

The above-mentioned test's number are calculated base on technical specifications of each analyzer. The real number of test per kit may higher than the calculation's number.

The above-mentioned test's number cover the loss of the dead volume of reagent bottles but not cover the loss of Calibrator and Control.

Please feel free to contact authorized distributor for further confirmation.

16. REFERENCES

- Adrian Bailey; Shamim S. Mohiuddin; Biochemistry, High Density Lipoprotein; 2020
- CLSI/NCCLS Evaluation of Precision Performance of Clinical Chemistry Devices, EP05-A2, 2004
- CLSI EPI7 · Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, 2nd Edition, 2017
- In house data, UMA Diagnostics

17. MANUFACTURER

UMA Co., Ltd.

2-19-6 Yokosuka, Matsudo City, Chiba

Prefecture 270-0031

TEL: 047-710-4871 (dial-in)

FAX: 047-710-4872