

**UMA CO., LTD.**

2-19-6 Yokosuka

Matsudo, Chiba, Japan



**MEASURE LDH IFCC**

Reagent for determination of Lactate Dehydrogenase

**IFCC Method**

↓ 2 - 8°C

IVD *In vitro* Diagnostics

QUALITY MANAGEMENT SYSTEM (BY TUV)

⊕ **DO NOT** freeze

⌚ 18 months/block from light

ISO 13485:2016

## 1. PURPOSE OF USE

Providing a quantitative *in vitro* assay for the Lactate Dehydrogenase (LDH) concentration in serum or plasma.

## 2. GENERAL INSTRUCTION

- For *in vitro* diagnostics use only.
- Diagnosis should be made in a comprehensive manner, in accordance with other related test results and clinical symptoms by the doctor in attendance.
- For guaranteed results, usage of this product must comply with the instruction in this manual.
- If you use automatic analyzers, follow their instructions carefully.

## SUMMARY

Lactate dehydrogenase (LDH) is an important enzyme of the anaerobic metabolic pathway. It belongs to the class of oxidoreductases, with an enzyme commission number EC 1.1.1.27. The function of the enzyme is to catalyze the reversible conversion of lactate to pyruvate with the reduction of NAD<sup>+</sup> to NADH and vice versa. The enzyme is present in a variety of organisms that include plants and animals. It is ubiquitously present in all tissues and serves as an important checkpoint of gluconeogenesis and DNA metabolism. A species-wide analysis of LDH demonstrates its well-preserved structure with only a few changes in the amino acid sequence across species. The structural similarity with slight amino acid changes provides a logical platform for designing functional molecules to modulate the catalytic potential and expression of the enzyme. In this article, we will focus our attention on the biochemical function, testing methods, and clinical relevance of the LDH enzyme.

## 3. MATERIALS REQUIRED BUT NOT INCLUDED

- Saline 0.9 % and high grade purified water
- Micropipet and other basic laboratory equipment.
- MEASURE Multi Calibrator and MEASURE Human Lyo L-1 and MEASURE Human Lyo L-2

## 4. REAGENT COMPOSITION & PREPARATION


- Reagent R-1: Lithium L-Lactate, Diethanolamine buffer  
Reagent R-1 is ready for use
- Reagent R-2: Nicotinamide Adenine Dinucleotide (NAD)  
Reagent R-2 is ready for use
- Once open, Reagent stored on board the instrument is stable for 30 days with Hitachi 7180 Analyzers.
- Applicable to various automated analyzers.
- Calibrator MEASURE Multi Calibrator (separately sold): Put 5 mL of purified water to the vials of Calibrator (MEASURE Multi Calibrator), leave at room temperature for 45 minutes and sometimes gently invert the vial before use. After reconstituting, Calibrator can be used without dilution.
- Controls MEASURE Human Lyo L-1 and MEASURE Human Lyo L-2 (separately sold): Put 5 mL of purified water to the vials of controls (Lyo L-1 and Lyo L-2); leave at room temperature for 45 minutes and sometimes gently invert the vial before use. After reconstituting, controls can be used without dilution.

## 5. SAMPLE PREPARATION & STORAGE

- Serum: Wait until the sample is completely coagulated. Take the supernatant to use as a specimen.
- Plasma: Treat blood sample by anticoagulant (Li-heparin and K2-EDTA); leave it to stand for 3 hours or centrifuge at 2000 rpm for 2 minutes; take the plasma layer (supernatant) and use as specimen.
- Analyze samples soon after collection. In case of storing sample 2 - 8°C, analyze within 7 days.
- Stability
  - 3 days at 15 - 25°C
  - 7 days at 2 - 8°C
  - 6 months at < -20°C
- See interferences section for details about possible sample interferences.


## 6. MEASUREMENT PRINCIPLE

Lactate Dehydrogenase (LDH) in sample is transformed into L-Lactate Pyruvic Acid under existence of Nicotinamide Adenine Dinucleotide (NAD). Obtain Lactate Dehydrogenase activity by measuring optically increase ratio of reduced form of Nicotinamide Adenine Dinucleotide (NADH) generated by this reaction.



## 7. ASSAY PROCEDURE

This product is compatible with various types of clinical analyzer. An example of the assay procedure is indicated below.



Perform the assay according to the instructions for operating the automated analyzer Hitachi models. Refer to the **13. INFORMATION FOR AUTOANALYZERS** for the details of the assay method. Contact **HUMA MEDICAL CO., LTD.** for information about the parameters for other automated analyzers.

## 8. CALCULATION & UNIT CONVERSION

### Calculation

- Calculate  $\Delta\text{Abs}$  of specimen & standards vs blank
- Plot a calibration curve  $\text{LDH} = f(\Delta\text{Abs})$
- Calculate LDH in specimen using the curve  
(doing same procedure for Controls)

### Unit conversion

$$\text{U/L} \times 0.0167 = \mu\text{kat/L}$$

## 9. PERFORMANCE & CORRELATION TEST

### a. Measuring range

- The assay is linear within a LDH enzyme activity of 3 - 1000 U/L.
- If the concentration of sample exceeds assay range, dilute the sample with saline and repeat the measurement.

### b. Detection Limit

Limit of Blank (LoB)	=	1.5 U/L
Limit of Detection (LoD)	=	3.0 U/L
Limit of Quantitation (LoQ)	=	3.0 U/L

The LoB, LoD and LoQ were determined in accordance with CLSI EP17-A2 requirements.

The LoB is the highest apparent analyte concentration expected to be found when replicates of a blank sample containing no analyte are tested. The LoB corresponds to the concentration below which analyte-free samples are found with a probability of 95%.

The LoD is determined based on the LoB and standard deviation of low concentration samples. The LoD corresponds to the lowest analyte concentration which can be detected (value above the LoB with a probability of 95%).

The LoQ is the lowest analyte concentration that can be reproducibly measured with a total error of 20%. It has been determined using low concentration samples.

### c. Performance

- Sensitivity: Change in absorbance when using purified water is below 0.01 Abs/min. and change in absorbance using solution (LDH 500 U/L) as sample is 0.03 - 0.07 Abs/minute.
- Accuracy: When measuring a control sample, the result is within  $\pm 10\%$  of assigned value.

### d. Precision (on Biolis 30i / SK300)

Representative performance data on the analyzers are given below.

Results obtained in individual laboratories may differ.

Precision was determined using controls followed the CLSI Approved Guideline EP5-A2 with repeatability, reproducibility and total precision (1 aliquot per run, 2 run per day, 20 days). The following results were obtained.

**Criterion: CV of Repeatability (aka. Within-run precision) is less than 3% and Total Precision is less than 5%.**



Measure LDH IFCC

Repeatability	Mean	SD	CV
	U/L	U/L	%
Control Lyo L-1	148.3	1.57	1.06
Control Lyo L-2	307.8	1.45	0.47

Reproducibility	Mean	SD	CV
	U/L	U/L	%
Control Lyo L-1	148.3	5.33	3.60
Control Lyo L-2	307.8	7.14	2.32

Total precision	Mean	SD	CV
	U/L	U/L	%
Control Lyo L-1	148.3	5.46	3.68
Control Lyo L-2	307.8	7.20	2.34

#### e. Correlation Test

Serum (n = 68)

Regression equation  $y = 1.01x + 1.50$

Correlation coefficient  $r = 0.999$

Plasma (n = 75)

Regression equation  $y = 1.01x - 1.61$

Correlation coefficient  $r = 1.000$

(y: value obtained from using UMA's reagent)

#### Reference Material for Calibration

ReCCS JCCLS CRM-001

### 10. EXPECTED VALUES

124 - 222 U/L

Reference range should be established at each facility and judgement should base on measurement results in a comprehensive manner together with clinical symptoms and other measurement results.

### 11. INTERFERENCES

- Icterus: No significant interference of conjugated/free bilirubin concentration up to 40 mg/dL

- Hemolysis: No significant interference of hemoglobin concentration in hemolyzed samples up to 18 mg/dL

- Lipemia (Intralipid): No significant interference triglycerides concentration up to 3000 FTU

- Ascorbic Acid: No significant interference of ascorbic acid concentration up to 50 mg/dL

- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings. Please use another methods if the result is affected by any factors

### 12. HANDLING, USAGE & DISPOSAL

#### Handling

1. Specimen can be potentially positive for infectious agents including hepatitis B virus and HIV. Wear glove and goggle when needed.

2. In case reagents got into skin, eye or mouth by mistake, wash it immediately with plenty of water and consult the doctor if needed.

3. If reagents are spilled, dilute with water and wipe it out. If specimen is spilled, spray 80% of alcohol over the specimen and wipe it out.

#### Usage

1. Store reagents under specified condition. Do not use after expiration date.

2. Do not use the container and auxiliaries included in this kit for other purposes.

3. Do not mix reagents of different lot for use.

4. Do not add to the reagent being used even if it is the same lot number.

#### Disposal

1. All specimens, as well as all instruments (e.g. test tubes) that come in contact with the specimens, must be treated by the following methods, or they must be treated according to the manual for infectious medical waste provided in each facility.

• Sterilize with an autoclave, subjecting them to high pressure saturated steam at 121 °C for more than 20 minutes. Do not process waste containing sodium hypochlorite solution with an autoclave.

• Immerse at least one hour in sodium hypochlorite solution (active chloride concentration of over 1000 ppm).

2. This reagent contains sodium azide. Sodium azide can react with lead pipe and/or steel pipe and can generate explosive metal azide. Make sure to use plenty of water at disposal. Concentration of sodium azide in R-2 is 0.05%.

**13. INFORMATION FOR AUTOANALYZERS**

## ❖ For Hitachi Model

Calculation Method		Rate
Temperature		37°C
	Specimen	4.0
Volume (μL)	R1	160
	R2	40
Wavelength (nm)	Main	340
	Sub	405
Measurement (cycle)	Point 1	10
	Point 2	21
	Point 3	28
Calibration type		Linear
Unit		U/L

**14. OTHER INSTRUCTIONS AND CAUTION**

- Results may differ depending on the sample/reagent ratio. Adjust parameters for different analyzer.
- Perform the QC procedure on the day of determination.

**15. PACKING AND KIT CONFIGURATION**

Code	Package	Test/Kit*	Test/Kit**
11L016A	1x60mL; 1x15mL	280	540
11L016A2	2x60mL; 2x15mL	560	1080
11L016A3	3x60mL; 3x15mL	840	1620
11L016A4	4x60mL; 4x15mL	1120	2160
11L006A	5x60mL; 5x15mL	1400	2700
11L016A6	6x60mL; 6x15mL	1680	3240
11L016	1x80mL; 1x20mL	380	750
11L016-2	2x80mL; 2x20mL	760	1500
11L006	3x80mL; 3x20mL	1140	2250
11L016-4	4x80mL; 4x20mL	1520	3000
11L016-5	5x80mL; 5x20mL	1900	3750

\* For middle-scale automatic analyzers such as: SK300; BS series; BA200; BA400. Chemwell Series; Dirui Series; Biolyzer series, HumanStar 300, Erba Series; Bioclab Series, BX 3010; Pictus P500,...

\*\* For large-scale automatic analyzers such as: CA800; CA400; Randox Imola; Randox Modena+; BM 6010; Biolis50i; SK500; AU Series; Pictus P700; C series; Ci series; HumanStar 600; Kenolab series ...

The above-mentioned test's number are calculated base on technical specifications of each analyzer. The real number of test per kit may higher than the calculation's number.

The above-mentioned test's number cover the loss of the dead volume of reagent bottles but not cover the loss of Calibrator and Control.

Please feel free to contact authorized distributor for further confirmation.

**16. REFERENCES**

1. Aisha Farhana; Sarah L. Lappin., Biochemistry, Lactate Dehydrogenase, 2020
1. CLSI/NCCLS Evaluation of Precision Performance of Clinical Chemistry Devices, EP05-A2, 2004
2. CLSI EP17 - Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, 2nd Edition, 2017
3. In house data, UMA Diagnostics

**17. MANUFACTURER**

UMA Co., Ltd.

2-19-6 Yokosuka, Matsudo City, Chiba

Prefecture 270-0031

TEL: 047-710-4871 (dial-in)

FAX: 047-710-4872