

UMA CO., LTD.

2-19-6 Yokosuka

Matsudo, Chiba, Japan



MEASURE TG(T)

Reagent for determination of Total Triglycerides

GK/GPO/POD Method

↓ 2 - 8°C

IVD *In vitro* Diagnostics

QUALITY MANAGEMENT SYSTEM (BY TUV)

⊛ DO NOT freeze

⌚ 18 months/block from light

ISO 13485:2016

1. PURPOSE OF USE

Providing a quantitative *in vitro* assay for the total triglycerides (TG(T)) concentration in serum or plasma.

2. GENERAL INSTRUCTION

- For *in vitro* diagnostics use only.
- Diagnosis should be made in a comprehensive manner, in accordance with other related test results and clinical symptoms by the doctor in attendance.
- For guaranteed results, usage of this product must comply with the instruction in this manual.
- If you use automatic analyzers, follow their instructions carefully.

SUMMARY

A triglyceride (TG) molecule consists of a glycerol backbone esterified with three fatty acids. Triglycerides are the main constituent of vegetable and animal fats in the diet, and are the main constituent of the body's fat stores. Serum or plasma total TG concentrations may be determined to assess metabolic disorders. TGs are not soluble in the aqueous environment of the bloodstream. For transport in the blood, they are carried by macromolecular particles called lipoproteins. The surface of lipoprotein particles is made up of proteins, free cholesterol, and phospholipids, oriented so as to be water soluble, and hydrophobic substances such as TG and esterified cholesterol are carried in the cores of the particles. Chylomicrons (CM) and very low density lipoprotein (VLDL) particles are the main TG carriers among the various types of lipoproteins, whereas low-density lipoproteins (LDL) and high-density lipoproteins (HDL) are mainly involved in transport of cholesterol. The CM particles carry dietary TG from the intestine to extrahepatic tissues such as muscle and adipose tissue. The VLDL particles carry TG made in the liver to extrahepatic tissues. Various disease states increase the amount of VLDL made by the liver and/or interfere with the clearance of CM or VLDL associated TG from the blood via uptake by muscle and/or adipose tissue. During conditions favoring lipolysis, nonesterified fatty acids (NEFAs) can be released from TG in the body's fat stores

and used for energy production. Excess mobilized NEFAs beyond what can be oxidized for energy may be reesterified into TG molecules in the liver. Accumulation of excess TG in the liver results in hepatic lipidosis (fatty liver). The equine liver is efficient at forming and secreting VLDL particles so as to remove excess TG, but disorders leading to impaired extrahepatic removal of the TG from those VLDL particles leads to hypertriglyceridemia, which, at levels greater than 400 to 500 mg/dL (4.52 - 5.65 mmol/L), causes lipemia.

3. MATERIALS REQUIRED BUT NOT INCLUDED

- Saline 0.9 % and high grade purified water
- Micropipet and other basic laboratory equipment.
- Lipids Calibrator and Lipids Control / MEASURE Human Lyo L-1 and MEASURE Human Lyo L-2

4. REAGENT COMPOSITION & PREPARATION

- Reagent R-1: LPL, POD, AOD, TOOS, N-ethyl-N-(2-hydroxy-3-sulfoethyl)-3-methylaniline sodium; Adenosine-5'-Triphosphate, Disodium (ATP 2Na); Glycerol Kinase (GK)

Reagent R-1 is ready for use

- Reagent R-2: GPO; Peroxidase (POD); 4-Aminoantipyrine (4-AA)

Reagent R-2 is ready for use

- Once open, Reagent stored on board the instrument is stable for 30 days with Hitachi 7180 Analyzers.

- Applicable to various automated analyzers.

- Calibrator Lipids Calibrator (separately sold): Put 2 mL of purified water to the vials of calibrator (Lipids Calibrator); leave at room temperature for 120 minutes and sometimes gently invert the vial before use. After reconstituting, Calibrator can be used without dilution.

- Control Lipids Control (separately sold): Put 2 mL of purified water to the vials of controls (Lipids Control); leave at room temperature for 120 minutes and sometimes gently invert the vial before use. After reconstituting, controls can be used without dilution.

Measure TG(T)

- Controls MEASURE Human Lyo L-1 and MEASURE Human Lyo L-2 (separately sold); Put 5 mL of purified water to the vials of controls (Lyo L-1 and Lyo L-2); leave at room temperature for 45 minutes and sometimes gently invert the vial before use. After reconstituting, controls can be used without dilution.

5. SAMPLE PREPARATION & STORAGE

- Serum: Wait until sample completely coagulated. Take the supernatant to use as specimen.

- Plasma: Collect blood after 12 - 14 hours fasting. Treat sample by anticoagulant (Li-heparin, K2-EDTA); leave sample to stand for 3 hours or centrifuge at 2000 rpm for 2 minutes; take the plasma layer (supernatant) and use as specimen.

- Analyze sample soon after collection.

- Stability


- 8 hours at 15 - 25°C
- 3 days at 2 - 8°C
- 30 days at < -20°C

- See interferences section for details about possible sample interferences.




6. MEASUREMENT PRINCIPLE

Free-glycerol is not eliminated in this Reagent. Measure the Total triglycerides. Ascorbate Oxidase in Reagent-1 eliminates reducing substance in the samples. Triglycerides in the samples are hydrolyzed by lipoprotein lipase in Reagent-2. In the next reaction, hydrogen peroxide is generated by glycerol kinase and glycerol-3-phosphate oxidase under the presence of ATP. This hydrogen peroxide induces oxidative condensation between 4-aminoantipyrin and TOOS, and generates red-violet quinone pigment. Concentration of triglycerides in the samples can be determined by measuring this absorbance.

1st reaction



2nd reaction



($\lambda_{\text{max}} = 555 \text{ nm}$)

7. ASSAY PROCEDURE

This product is compatible with various types of clinical analyzer. An example of the assay procedure is indicated below.



————→ ΔAbs at 546/800nm ———→ Concentration TG(T)

Perform the assay according to the instructions for operating the automated analyzer Hitachi models. Refer to the **12. INFORMATION FOR AUTOANALYZERS** for the details of the assay method. Contact **HUMA MEDICAL CO., LTD.** for information about the parameters for other automated analyzers.

8. CALCULATION & UNIT CONVERSION

Calculation

- Calculate ΔAbs of specimen & standards vs blank

- Plot a calibration curve $\text{TG(T)} = f(\Delta\text{Abs})$

- Calculate TG(T) in specimen using the curve

(doing same procedure for Controls)

Unit conversion

$$\text{mg/dL} \times 0.0113 = \text{mmol/L}$$

9. PERFORMANCE & CORRELATION TEST

a. Measuring range

- The assay is linear within an TG(T) concentration range in serum/plasma of 0.11 - 11.3 mmol/L.

- If the concentration of sample exceeds assay range, dilute the sample with saline and repeat the measurement.

b. Detection Limit

Limit of Blank (LoB)	=	0.01 mmol/L
Limit of Detection (LoD)	=	0.03 mmol/L
Limit of Quantitation (LoQ)	=	0.11 mmol/L

The LoB, LoD and LoQ were determined in accordance with CLSI EP17-A2 requirements.

The LoB is the highest apparent analyte concentration expected to be found when replicates of a blank sample containing no analyte are tested. The LoB corresponds to the concentration below which analyte-free samples are found with a probability of 95%.

The LoD is determined based on the LoB and standard deviation of low concentration samples. The LoD corresponds to the lowest analyte concentration which can be detected (value above the LoB with a probability of 95%).

The LoQ is the lowest analyte concentration that can be reproducibly measured with a total error of 20%. It has been determined using low concentration samples.

c. Performance

- Sensitivity: When using purified water, absorbance change is 0.001 - 0.050; using triglyceride 1.13 mmol/L solution, absorbance change is 0.10 - 0.50.

- Accuracy: When measuring a control sample, the result is within $\pm 10\%$ of assigned value.

d. Precision (on Biolis 30i / SK300)

Representative performance data on the analyzers are given below.

Results obtained in individual laboratories may differ.

Precision was determined using controls followed the CLSI Approved Guideline EP5-A2 with repeatability, reproducibility and total precision (1 aliquot per run, 2 run per day, 20 days). The following results were obtained.

Criterion: CV of Repeatability (aka. Within-run precision) is less than 3% and Total Precision is less than 5%.

Repeatability	Mean mmol/L	SD mmol/L	CV %
Control Lyo L-1	1.08	0.02	1.65
Control Lyo L-2	4.60	0.04	0.93

Reproducibility	Mean mmol/L	SD mmol/L	CV %
Control Lyo L-1	1.08	0.03	2.50
Control Lyo L-2	4.60	0.07	1.48

Total precision	Mean mmol/L	SD mmol/L	CV %
Control Lyo L-1	1.08	0.03	2.76
Control Lyo L-2	4.60	0.07	1.62

e. Correlation Test

Same measuring principle

Serum

Regression equation: $y = 0.9692x - 2.50$ (n = 50)

Correlation coefficient: $r = 0.9992$

Plasma

Regression equation: $y = 1.0210x - 1.00$ (n = 50)

Correlation coefficient: $r = 0.9993$

(y: value obtained from using UMA's reagent)

Reference Materials for Calibration

ReCCS JCCRM 223

10. EXPECTED VALUES

Normal reference range

- Male 0.45 - 2.65 mmol/L

- Female 0.34 - 1.33 mmol/L

Reference range should be established at each facility and judgement should be based on measurement results in a comprehensive manner together with clinical symptoms and other measurement results.

11. INTERFERENCES

- Icterus: No significant interference of conjugated/free bilirubin concentration up to 20 mg/dL

- Hemolysis: No significant interference of hemoglobin concentration up to 200 mg/dL

- Ascorbic Acid: No significant interference of ascorbic acid concentration up to 50 mg/dL

- For diagnostic purposes, the results should always be assessed in conjunction with the patient's medical history, clinical examination and other findings. Please use another methods if the result is affected by any factors

12. INFORMATION FOR AUTOANALYZERS

❖ For Hitachi Model		
Calculation Method	Two point	
Temperature	37°C	
Specimen	2.0	
Volume (µL)	R1	150
	R2	50
Wavelength (nm)	Main	546
	Sub	800
Measurement (cycle)	Point 1	10
	Point 2	16
	Point 3	34

Calibration type	Linear
Unit	mg/dL

13. HANDLING, USAGE & DISPOSAL

Handling

1. Specimen can be potentially positive for infectious agents including hepatitis B virus and HIV. Wear glove and goggle when needed.

2. In case reagents got into skin, eye or mouth by mistake, wash it immediately with plenty of water and consult the doctor if needed.

3. If reagents are spilled, dilute with water and wipe it out. If specimen is spilled, spray 80% of alcohol over the specimen and wipe it out.

Usage

1. Store reagents under specified condition. Do not use after expiration date.

2. Do not use the container and auxiliaries included in this kit for other purposes.

3. Do not mix reagents of different lot for use.

4. Do not add to the reagent being used even if it is the same lot number.

Disposal

1. All specimens, as well as all instruments (e.g. test tubes) that come in contact with the specimens, must be treated by the following methods, or they must be treated according to the manual for infectious medical waste provided in each facility.

- Sterilize with an autoclave, subjecting them to high pressure saturated steam at 121 °C for more than 20 minutes. Do not process waste containing sodium hypochlorite solution with an autoclave.

- Immerse at least one hour in sodium hypochlorite solution (active chloride concentration of over 1000 ppm).

2. This reagent contains sodium azide. Sodium azide can react with lead pipe and/or steel pipe and can generate explosive metal azide. Make sure to use plenty of water at disposal. Concentration of sodium azide in R-2 is 0.05%.

14. OTHER INSTRUCTIONS AND CAUTION

- Results may differ depending on the sample/reagent ratio. Adjust parameters for different analyzer.

- Perform the QC procedure on the day of determination.

15. PACKING AND KIT CONFIGURATION

Code	Package	Test/Kit*	Test/Kit**
11T015A	1x60mL; 1x20mL	310	540

11T015A2	2x60mL; 2x20mL	620	1080
11T015A3	3x60mL; 3x20mL	930	1620
11T015A4	4x60mL; 4x20mL	1240	2160
11T005A	5x60mL; 5x20mL	1550	2700
11T015A6	6x60mL; 6x20mL	1860	3240
11T015	1x90mL; 1x30mL	470	810
11T015-2	2x90mL; 2x30mL	940	1620
11T005	3x90mL; 3x30mL	1410	2430
11T015-4	4x90mL; 4x30mL	1880	3240
11T015-5	5x90mL; 5x30mL	2350	4050

* For middle-scale automatic analyzers such as: SK300; BS series; BA200; BA400. Chemwell Series; Dirui Series; Biolyzer series, HumanStar 300, Erba Series; Bioelab Series, BX 3010; Pictus P500;...

** For large-scale automatic analyzers such as: CA800; CA400; Randox Imola; Randox Modena+; BM 6010; Biolis50i; SK500; AU Series; Pictus P700; C series; Ci series; HumanStar 600; Kenolab series ...

The above-mentioned test's number are calculated base on technical specifications of each analyzer. The real number of test per kit may higher than the calculation's number.

The above-mentioned test's number cover the loss of the dead volume of reagent bottles but not cover the loss of Calibrator and Control.

Please feel free to contact authorized distributor for further confirmation.

16. REFERENCES

1. MARLYN S. WHITNEY, In Clinical Veterinary Advisor: The Horse, 2012
2. CLSI/NCCLS Evaluation of Precision Performance of Clinical Chemistry Devices, EP05-A2, 2004
3. CLSI EP17 · Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures, 2nd Edition, 2017
4. In house data, UMA Diagnostics

17. MANUFACTURER

UMA Co., Ltd.

2-19-6 Yokosuka, Matsudo City, Chiba

Prefecture 270-0031

TEL: 047-710-4871 (dial-in)

FAX: 047-710-4872