

試験成績書
Certificate of Analysis

株式会社ユーエムエー
UMA CO.,LTD

試験書番号 Certificate No.	1067
コード Code No.	11P002
製品名 Product	メジャー PG II (3:1)
試験実施年月日 Examination date	Year MM DD 2023/2/8
総合判定 Result	合格

Examination Reagent	R-1:Code No. 11P102	Lot No.	164ACB
	R-2:Code No. 11P202	Lot No.	164ACB

Manufactured Date	2023/2/8
	Year MM DD

Expiry Date 有効期限	2024/08
	Year MM DD

Analytes			Standard		Result
Appearance	R-1		Colorless		Pass test
	R-2		Opalescent		Pass test
pH (at 20℃)	R-1		7.00	~ 8.00	7.4
	R-2		6.00	~ 8.00	6.2
Blank		mAbs.	less than 20		-2
Sensitivity	Cal.Level-6	mAbs.	more than 500		655
Precision /with in run (C.V.)			≤ 10.0% (n=10)		0.48%
* Accuracy (ng/mL)	Control-1	19.5	90%	~ 110%	100%
	Control-2	68.4	90%	~ 110%	97%

*コントロールのロット変更時は値を変更する。

試験者	品質管理課長	製造管理者