

試験成績書
Certificate of Analysis

株式会社ユーエムエー
UMA CO.,LTD

試験書番号 Certificate No.	1070
コード Code No.	11P002
製品名 Product	メジャー PG II (3:1)
試験実施年月日 Examination date	Year MM DD 2024/3/7
総合判定 Result	合格

Examination Reagent	R-1:Code No. 11P102	Lot No.	168ADC
	R-2:Code No. 11P202	Lot No.	168ADC

Manufactured Date	2024/3/7
	Year MM DD

Expiry Date 有効期限	2025/09
	Year MM DD

Analytes			Standard			Result
Appearance	R-1		Colorless			Pass test
	R-2		Opalescent			Pass test
pH (at 20℃)	R-1		7.00	~	8.00	7.4
	R-2		6.00	~	8.00	6.2
Blank		mAbs.	less than 20			-4
Sensitivity	Cal.Level-6	mAbs.	more than 500			576
Precision /with in run (C.V.)			≤ 10.0% (n=10)			0.95%
*Accuracy (ng/mL)	Control-1	19.2	90%	~	110%	103%
	Control-2	66.9	90%	~	110%	102%

*コントロールのロット変更時は値を変更する。

試験者	品質管理課長	製造管理者
		