

試験成績書
Certificate of Analysis

株式会社ユーエムエー
UMA CO.,LTD

試験書番号 Certificate No.	1066
コード Code No.	11P001
製品名 Product	メジャー PG I (3:1)
試験実施年月日 Examination date	Year MM DD 2023/2/8
総合判定 Result	合格

Examination Reagent	R-1:Code No. 11P101	Lot No.	164ACB
	R-2:Code No. 11P201	Lot No.	164ACB

Manufactured Date

2023/2/8
Year MM DD

Expiry Date
有効期限

2024/08
Year MM DD

Analytes		Standard		Result
Appearance	R-1	Colorless		Pass test
	R-2	Opalescent		Pass test
pH (at 20℃)	R-1	7.00 ~ 8.00		7.4
	R-2	6.00 ~ 8.00		7.4
Blank	mAbs.	less than 20		-4
Sensitivity	Cal.Level-6 mAbs.	more than 500		643
Precision /with in run (C.V.)		≤ 10.0% (n=10)		0.81%
* Accuracy (ng/mL)	Control-1 49.2	90% ~ 110%		97%
	Control-2 156.3	90% ~ 110%		102%

*コントロールのロット変更時は値を変更する。

試験者	品質管理課長	製造管理者
		