

Bilirubin Calibrator Certificate of Analysis

Issued Date : 2023/5/1
 Certificate No. : 123
 Lot No. : 009ACA
 Manufacture date : 2023/5/1
 Expiration date : 2024/11


Origin : Freeze Dried Pooled Human Serum
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for bilirubin Calibrator

Analyte	Standard	Result
T-BIL		12.0 mg/dL
D-BIL		9.4 mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
HBV Ag.Ab	Negative (—)	Negative (—)
HCV Ag.Ab	Negative (—)	Negative (—)
HIV Ag.Ab	Negative (—)	Negative (—)
TP	Negative (—)	Negative (—)
ATLV Ab.	Negative (—)	Negative (—)

Preparation : Reconstitute 2 mL of distilled water and
 then stand for 20 minutes, swirling occasionally.
 Store in the dark after reconstitution.

Measurement by :
 Date :


 2023/5/1