

# Bilirubin Calibrator Certificate of Analysis

Issued Date : 2024/3/1  
Certificate No. : 128  
Lot No. : 009ADC  
Expiration date : 2025/09

Origin : Freeze Dried Pooled Human Serum  
each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for bilirubin Calibrator

Analyte	Standard	Result
T-BIL		12.0 mg/dL
D-BIL		9.4 mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
HBV Ag.Ab	Negative (－)	Negative (－)
HCV Ag.Ab	Negative (－)	Negative (－)
HIV Ag.Ab	Negative (－)	Negative (－)
Tp	Negative (－)	Negative (－)
ATLV Ab.	Negative (－)	Negative (－)

Preparation : Reconstitute 2 mL of distilled water and then stand for 20 minutes, swirling occasionally. Store in the dark after reconstitution.

Measurement by : Takeshi Terahara  
Date : 2024/3/1