

Lipids Calibrator Certificate of Analysis

Issued Date : 2024/3/4
Certificate No. : 317
Lot No. : 120ADC
Expiration date : 2025/09

Origin : Freeze Dried Pooled Human Serum
each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte		Standard		Result	
HDL-C	HDL	≥50 mg/dL	59.0 mg/dL	1.53	mmol/L
LDL-C	LDL		136 mg/dL	3.52	mmol/L
TC			210 mg/dL	5.43	mmol/L
TG			118 mg/dL	1.33	mmol/L
T-TG			120 mg/dL	1.35	mmol/L
HBV Ag.Ab		Negative (－)	Negative (－)		
HCV Ag.Ab		Negative (－)	Negative (－)		
HIV Ag.Ab		Negative (－)	Negative (－)		
Tp		Negative (－)	Negative (－)		
ATLV Ab.		Negative (－)	Negative (－)		

Preparation : Reconstitute 2 mL of distilled water and mix gently,
then store room temperature and sometimes mix
during 2 hours before use.

Measurement by : Takeshi Yoshikawa
Date : 2024/3/4