

Lipids Calibrator Certificate of Analysis

Issued Date : 2024/5/1
Certificate No. : 320
Lot No. : 120ADE
Expiration date : 2025/11

Origin : Freeze Dried Pooled Human Serum
each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

| Analyte | | Standard | | Result | |
|-----------|-----|--------------|-------|--------------|--------|
| HDL-C | HDL | ≥50 | mg/dL | 55.0 | mg/dL |
| | | | | 1.42 | mmol/L |
| LDL-C | LDL | | | 113 | mg/dL |
| | | | | 2.92 | mmol/L |
| TC | | | | 208 | mg/dL |
| | | | | 5.38 | mmol/L |
| TG | | | | 131 | mg/dL |
| T-TG | | | | 131 | mg/dL |
| | | | | 1.48 | mmol/L |
| | | | | 1.48 | mmol/L |
| HBV Ag.Ab | | Negative (－) | | Negative (－) | |
| HCV Ag.Ab | | Negative (－) | | Negative (－) | |
| HIV Ag.Ab | | Negative (－) | | Negative (－) | |
| Tp | | Negative (－) | | Negative (－) | |
| ATLV Ab. | | Negative (－) | | Negative (－) | |

Preparation : Reconstitute 2 mL of distilled water and mix gently,
then store room temperature and sometimes mix
during 2 hours before use.

Measurement by : Takeshi Yoshino
Date : 2024/5/1