

PG I /PG II Control

PG I /PG II Control Set

Intended Use

Control for PG I /PG II

Precautions

HBs Ag :	Negative	(-)
HCV Ab :	Negative	(-)
HIV Ab :	Negative	(-)

The product must however be handled as potentially biohazardous material.

This product includes 0.09w/v% sodium azide as preservative.

After the measurement, it must be diluted with plenty of water and then drained away.

Storage : Dark and Cold (2~8 °C)

Lot No. : 127ABI

Manufacture Date : 2022/9/29

Exp. Date : 2024.03.

【Control value】

Level	Color (cap)	Measured value (ng/mL)	Reference range
PG I	LOW Purple	50.7	45.7 - 55.7
	High Brown	152.8	137.8 - 167.8
PG II	LOW Purple	19.6	17.6 - 21.6
	High Brown	70.5	63.5 - 77.5

manufactured by :

UMA Co., Ltd.

2-19-6 Yokosuka Matsudo City Chiba-Ken,

270-0031, Japan

phone: +81-47-710-4871

試験者	品質管理課長	製造管理者
		