

## PG I /PG II Control

## PG I /PG II Control Set

## Intended Use

Control for PG I /PG II

## Precautions

HBs Ag :	Negative	(-)
HCV Ab :	Negative	(-)
HIV Ab :	Negative	(-)

The product must however be handled as potentially biohazardous material.

This product includes 0.09w/v% sodium azide as preservative.

After the measurement, it must be diluted with plenty of water and then drained away.

Storage : Dark and Cold (2~8 °C)

Lot No. : 131ACL

Exp.Date : 2025.06.

## 【Control value】

Level	Color (cap)	Measured value (ng/mL)	Reference range
PG I	LOW Purple	47.8	42.8 - 52.8
	High Brown	145.3	135.3 - 155.3
PG II	LOW Purple	18.6	16.6 - 20.6
	High Brown	65.6	58.6 - 72.6

manufactured by :

UMA Co., Ltd.

2-19-6 Yokosuka Matsudo City Chiba-Ken,

270-0031, Japan

phone: +81-47-710-4871

試験者	品質管理課長	製造管理者
		