

## Lipids Calibrator Certificate of Analysis

Issued Date : 2024/8/1  
Certificate No. : 326  
Lot No. : 121ADH  
Expiration date : 2026/02

Origin : Freeze Dried Pooled Human Serum  
each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte		Standard		Result			
HDL-C	HDL	$\geq 50$	mg/dL	55.0	mg/dL	1.42	mmol/L
LDL-C	LDL			113	mg/dL	2.92	mmol/L
TC				208	mg/dL	5.38	mmol/L
TG				131	mg/dL	1.48	mmol/L
T-TG				131	mg/dL	1.48	mmol/L
HBV Ag.Ab		Negative (—)		Negative (—)			
HCV Ag.Ab		Negative (—)		Negative (—)			
HIV Ag.Ab		Negative (—)		Negative (—)			
Tp		Negative (—)		Negative (—)			
ATLV Ab.		Negative (—)		Negative (—)			

Preparation : Reconstitute 2 mL of distilled water and mix gently,  
then store room temperature and sometimes mix  
during 2 hours before use.

Measurement by : Takeshi Yoshitake  
Date : 2024/8/1