

# Lipids Calibrator Certificate of Analysis

Issued Date : 2024/10/1  
 Certificate No. : 328  
 Lot No. : 121ADJ  
 Expiration date : 2026/04

Origin : Freeze Dried Pooled Human Serum  
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte		Standard		Result	
HDL-C	HDL	$\geq 50$ mg/dL	55.0 mg/dL	1.42	mmol/L
LDL-C	LDL		113 mg/dL	2.92	mmol/L
TC			208 mg/dL	5.38	mmol/L
TG			131 mg/dL	1.48	mmol/L
T-TG			131 mg/dL	1.48	mmol/L
HBV Ag.Ab		Negative (—)	Negative (—)		
HCV Ag.Ab		Negative (—)	Negative (—)		
HIV Ag.Ab		Negative (—)	Negative (—)		
Tp		Negative (—)	Negative (—)		
ATLV Ab.		Negative (—)	Negative (—)		

**Preparation** : Reconstitute 2 mL of distilled water and mix gently,  
 then store room temperature and sometimes mix  
 during 2 hours before use.

Measurement by : Takekazu Yukioka  
 Date : 2024/10/1