

# Lipids Calibrator Certificate of Analysis

Issued Date : 2021/8/1  
 Certificate No. : 284  
 Lot No. : 117AAH  
 Expiration date : 2023/02

Origin : Freeze Dried Pooled Human Serum  
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte		Standard	Result			
HDL-C	HDL	$\geq 50$ mg/dL	66.7	mg/dL	1.73	mmol/L
LDL-C	LDL		135	mg/dL	3.50	mmol/L
TC			213	mg/dL	5.52	mmol/L
TG			118	mg/dL	1.33	mmol/L
T-TG			125	mg/dL	1.41	mmol/L
HBV Ag.Ab		Negative (—)	Negative (—)			
HCV Ag.Ab		Negative (—)	Negative (—)			
HIV Ag.Ab		Negative (—)	Negative (—)			
Tp		Negative (—)	Negative (—)			
ATLV Ab.		Negative (—)	Negative (—)			

**Preparation** : Reconstitute 2 mL of distilled water and mix gently,  
 then store room temperature and sometimes mix  
 during 2 hours before use.

Measurement by : Takeshi Tsubota  
 Date : 2021/8/1