

hs-CRP Calibrator Set

use for Measure hs-CRP reagent

Intended Use

This product is for purpose of calibrating results in the quantitative determination of C-reactive protein in Serum or Plasma automated immunoassay. For in vitro diagnostic use only.

Precautions

HBs Ag : Negative (—)
HCV Ab : Negative (—)
HIV Ab : Negative (—)

The product must however be handled as potentially biohazardous material.

This product includes 0.09w/v% sodium azide as preservative.

In case the product gets into eye or mouth, or comes into contact with skin, wash with plenty of running water and consult a doctor if needed.

Storage : Dark and Cold (2~8 °C)

Lot No. : 120ABJ

Exp. Date : 2024. 04.

【Assigned value】

| Level | Color (cap) | mg/dL | mg/L |
|-------|-------------|-------|------|
| 1 | Colorless | 0.00 | 0.0 |
| 2 | Blue | 0.16 | 1.6 |
| 3 | Green | 1.65 | 16.5 |
| 4 | Yellow | 4.40 | 44.0 |
| 5 | Orange | 22.7 | 227 |
| 6 | Pink | 42.4 | 424 |

Manufactured by :

UMA Co., Ltd.

2-19-6 Yokosuka Matsudo City Chiba-Ken,

270-0031 Japan

phone: +81-47-710-4871

| 試験者 | 品質管理課長 | 製造管理者 |
|---|---|---|
|  |  |  |