

## Lipids Calibrator Certificate of Analysis

Issued Date : 2024/3/4  
 Certificate No. : 317  
 Lot No. : 120ADC  
 Expiration date : 2025/09

Origin : Freeze Dried Pooled Human Serum  
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte		Standard	Result	
HDL-C	HDL	$\geq 50$ mg/dL	59.0 mg/dL	1.53 mmol/L
LDL-C	LDL		136 mg/dL	3.52 mmol/L
TC			210 mg/dL	5.43 mmol/L
TG			118 mg/dL	1.33 mmol/L
T-TG			120 mg/dL	1.35 mmol/L
HBV Ag.Ab		Negative (-)	Negative (-)	
HCV Ag.Ab		Negative (-)	Negative (-)	
HIV Ag.Ab		Negative (-)	Negative (-)	
TP		Negative (-)	Negative (-)	
ATLV Ab.		Negative (-)	Negative (-)	

Preparation : Reconstitute 2 mL of distilled water and mix gently,  
 then store room temperature and sometimes mix  
 during 2 hours before use.

Measurement by :  
 Date :

*Takekazu Yoshino*  
 2024/3/4