

Lipids Control Certificate of Analysis

Issued Date : 2024/12/2
 Certificate No. : 332
 Lot No. : 122ADL
 Expiration date : 2026/06

Origin : Freeze Dried Pooled Human Serum
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte	Standard	Result	(Reference Range)
HDL-C	HDL	52.0	46.0-58.0 mg/dL
		1.34	1.19-1.50 mmol/L
LDL-C	LDL	111.0	98.0-124 mg/dL
		2.87	2.53-3.21 mmol/L
TC		196.0	174-218 mg/dL
		5.07	4.50-5.64 mmol/L
TG		118.0	106-130 mg/dL
		1.33	1.20-1.47 mmol/L
(T-TG)		(121)	(109-133) mg/dL
		1.37	(1.23-1.50) mmol/L
HBV Ag.Ab	Negative (-)	Negative (-)	
HCV Ag.Ab	Negative (-)	Negative (-)	
HIV Ag.Ab	Negative (-)	Negative (-)	
Tp	Negative (-)	Negative (-)	
ATLV Ab.	Negative (-)	Negative (-)	

Preparation : Reconstitute 2 mL of distilled water and mix gently, then store room temperature and sometimes mix during 2 hours before use.

Measurement by :
 Date :

Takeshi Yokokawa

 2024/12/2
