

Bilirubin Calibrator Certificate of Analysis

Issued Date : 2025/6/2
 Certificate No. : 14
 Lot No. : 009AEF
 Expiration date : 2026/12

Origin : Freeze Dried Pooled Human Serum
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for bilirubin Calibrator

Analyte	Standard	Result
T-BIL		12.0 mg/dL
D-BIL		9.4 mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
		mg/dL
HBV Ag.Ab	Negative (-)	Negative (-)
HCV Ag.Ab	Negative (-)	Negative (-)
HIV Ag.Ab	Negative (-)	Negative (-)
Tp	Negative (-)	Negative (-)
ATLV Ab.	Negative (-)	Negative (-)

Preparation : Reconstitute 2 mL of distilled water and then stand for 20 minutes, swirling occasionally. Store in the dark after reconstitution.

Measurement by : Takeshi Yoshida
 Date : 2025/6/2