

# Bilirubin Calibrator Certificate of Analysis

Issued Date : 2026/2/2  
 Certificate No. : 16  
 Lot No. : 009AFB  
 Manufacture date : 2026/2/2  
 Expiration date : 2027/08

Origin : Freeze Dried Pooled Human Serum  
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for bilirubin Calibrator

| Analyte   | Standard     | Result       |
|-----------|--------------|--------------|
| T-BIL     |              | 12.0 mg/dL   |
| D-BIL     |              | 9.4 mg/dL    |
|           |              | mg/dL        |
|           |              | mg/dL        |
|           |              | mg/dL        |
|           |              | mg/dL        |
|           |              | mg/dL        |
|           |              | mg/dL        |
| HBV Ag.Ab | Negative (-) | Negative (-) |
| HCV Ag.Ab | Negative (-) | Negative (-) |
| HIV Ag.Ab | Negative (-) | Negative (-) |
| Tp        | Negative (-) | Negative (-) |
| ATLV Ab.  | Negative (-) | Negative (-) |

Preparation : Reconstitute 2 mL of distilled water and then stand for 20 minutes, swirling occasionally. Store in the dark after reconstitution.

Measurement by : Takeshi Yuhiko  
 Date : 2026/2/2