

Lipid Calibrator Certificate of Analysis

Issued Date : 2025/12/1
 Certificate No. : 342
 Lot No. : 123AEL
 Expiration date : 2027/06

Origin : Freeze Dried Pooled Human Serum
 each specimen has passed HBV , HCV and HIV tests.

Table 1, Certified Concentrations for Lipid Calibrator

Analyte		Standard	Result	
HDL-C	HDL	≥ 50 mg/dL	64.0 mg/dL	1.66 mmol/L
LDL-C	LDL		116 mg/dL	3.00 mmol/L
TC			210 mg/dL	5.43 mmol/L
TG			124 mg/dL	1.40 mmol/L
T-TG			125 mg/dL	1.41 mmol/L
HBV Ag.Ab		Negative (-)	Negative (-)	
HCV Ag.Ab		Negative (-)	Negative (-)	
HIV Ag.Ab		Negative (-)	Negative (-)	
TP		Negative (-)	Negative (-)	
ATLV Ab.		Negative (-)	Negative (-)	

Preparation : Reconstitute 2 mL of distilled water and mix gently,
 then store room temperature and sometimes mix
 during 2 hours before use.

Measurement by : Takashi Yoshida
 Date : 2025/12/1

UMA Co., Ltd. :2-19-6 Yokosuka Matsudo City Chiba-Ken,270-0031, Japan